


# Stealth Cage

The BritzCo Stealth Cage is the perfect addition to Blackrock's CerePlex W wireless system. Benefit from optimal wireless transmission regardless of where the antennae are exactly placed. CerePlex W and the BritzCo Stealth Cage are the integrated solution for interference free wireless neural recordings. The cage has two defined cells each with multiple doors for ease of access to any part of the inside of the cage.

## Perfect match for CerePlex W wireless recording system

### Features:

- 95% of this cage is formed with non-metallic materials.
- Cage case is fabricated from high pressure compact laminates that are known to be x-ray and RF transparent.
- Polycarbonates are used for clear panels, while polyethylene is used for sliding surfaces.
- Front panels are formed with protruded fiberglass and thermo set polyester bars and rods.



### Specifications:

- 73" tall
- 32" deep
- 31" wide
- Each unit is about 31" tall

For more information, please contact us: [sales@blackrockmicro.com](mailto:sales@blackrockmicro.com)